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Referenzen

1. Korniewicz DM, Kirwin M, Cresci K, Larson E. Leakage of latex and vinyl exam gloves in high and low risk clinical settings. *Am Ind Hyg Assoc J* 1993;54(1):22-26
2. Douglas A, Simon TR, Goddard M. Barrier durability of latex and vinyl medical gloves in clinical settings. *Am Ind Hyg Assoc J* 1997;58(9):672-676
3. Rego A, Roley L. In-use barrier integrity of gloves: latex and nitrile superior to vinyl. *Am J Infect Control* 1999;27(5):405-410
4. Korniewicz DM, El-Masri M, Broyles JM, Martin CD, O'connell KP. Performance of latex and nonlatex medical examination gloves during simulated use. *Am J Infect Control* 2002;30(2):133-138
5. Kerr LN, Chaput MP, Cash LD, O'Malley LG, Sarhrani EM, Teixeira JC, Boivin WS, Mailhot SA. Assessment of the durability of medical examination gloves. *J Occup Environ Hyg* 2004;1(9):607-612
6. Korniewicz DM, Kirwin M, Cresci K, Sing T, Choo TE, Wool M, Larson E. Barrier protection with examination gloves: double versus single. *Am J Infect Control* 1994;22(1):12-15
7. Korniewicz DM, Laughon BE, Butz A, Larson E. Integrity of vinyl and latex procedure gloves. *Nurs Res* 1989;38(3):144-146
8. Korniewicz DM, Laughon BE, Cyr WH, Lytle CD, Larson E. Leakage of virus through used vinyl and latex examination gloves. *J Clin Microbiol* 1990;28(4):787-788
9. Olsen RJ, Lynch P, Coyle MB, Cummings J, Bokete T, Stamm WE. Examination gloves as barriers to hand contamination in clinical practice. *JAMA* 1993;270(3):350-353
10. Gerhardt GG. Results of microbiological investigations on the permeability of procedure and surgical gloves. *Zentralbl Hyg Umwelt-med* 1989;188(3-4):336-342
11. Klein RC, Party E, Gershey EL. Virus penetration of examination gloves. *Biotechniques* 1990;9(2):196-199
12. Neal JG, Jackson EM, Suber F, Edlich R. Latex glove penetration by pathogens: a review of the literature. *J Long Term Eff Med Implants* 1998;8(3-4):233-240
13. AAMI. Safe use and handling of glutaraldehyde based products in health care facilities. American National Standard 1996
14. Boyce JM, Pittet D. Guideline for hand hygiene in health-care settings. Recommendations of the healthcare infection control practices advisory committee and the HICPAC/SHEA/APIC/IDSA hand hygiene task force. *MMWR Recomm Rep* 2002;51:1-45
15. Laidlaw JL, Connor TH, Theiss JC, Anderson RW, Matney TS. Permeability of latex and polyvinyl chloride gloves to 20 antineoplastic drugs. *Am J Hosp Pharm* 1984;41(12):2618-2623
16. Walleamacq PE, Capron A, Vanbinst R, Boeckmans E, Gillard J, Favier B. Permeability of 13 different gloves to 13 cytotoxic agents under controlled dynamic conditions. *Am J Health Syst Pharm* 2006;63(6):547-556
17. Johnson F. Disposable gloves: research findings on use in practice. *Nurs Stand* 1997;11(16):39-40
18. Burke FJ, Watts DC, Wilson NH. Some physical factors influencing tactile perception with disposable non-sterile gloves. *J Dent* 1989;17(2):72-76
19. Infection Control Nurse Association (ICNA). *Gloves Uses Guidelines*, UK, September 1999
20. Hunte SC. Choosing the right glove for the right purpose. *Prof Nurse* 2004;20(3):43-47
21. Siegel JD, Rhinehart E, Jackson M, Chiarello L, Healthcare Infection Control Practice Advisory Committee. Guideline for isolation precautions: preventing transmission of infectious agents in healthcare settings. *Am J Infect Control* 2007;35(Suppl 2):S65-S164
22. Aalto-Korte K, Alanko K, Henriks-Eckerman ML, Estlander T, Jolanki R. Allergic contact dermatitis from bisphenol A in PVC gloves. *Contact Dermatitis* 2003;49(4):202-205
23. Sowa J, Kobayashi H, Tsuruta D, Sugawara K, Ishii M. Allergic contact dermatitis due to adipic polyester in vinyl chloride gloves. *Contact Dermatitis* 2005;53(4):243-244
24. Aalto-Korte K, Alanko K, Henriks-Eckerman ML, Jolanki R. Antimicrobial allergy from polyvinyl chloride gloves. *Arch Dermatol* 2006;142(10):1326-1330
25. Aalto-Korte K, Ackermann L, Henriks-Eckerman ML, Välimaa J, Reinikka-Railo H, Leppänen E, Jolanki R. 1,2-Benzisothiazolin-3-one in disposable polyvinyl chloride gloves for medical use. *Contact Dermatitis* 2007;57(6):365-370
26. Matthieu L, Godoi AF, Lambert J, Van Grieken R. Occupational allergic contact dermatitis from bisphenol A in vinyl gloves. *Contact Dermatitis* 2003; 49(6):281-283
27. Ueno M, Adachi A, Horikawa T, Inoue N, Mori A, Sasaki K. Allergic contact dermatitis caused by poly(adipic acid-co-1,2-propylene glycol) and di-(n-octyl) tin-bis (2-ethylhexyl maleate) in vinyl chloride gloves. *Contact Dermatitis* 2007;57(5):349-351
28. Ban of phthalates in childcare articles and toys, press release IP/99/829, 10 November 1999; http://europa.eu/index_en.htm IP/99/829
29. Brown P, KrennHrubec K. Phthalates and Children's Products; <http://www.center4research.org/phthalates.html>, July 2008
30. Substances of Very High Concern: Annex XV reports to be commented by Interested Parties; http://echa.europa.eu/consultations/authorisation/svhc/svhc_cons_en.asp
31. European Directive 2007/19/EC, Official Journal of the European Union L 97. 12 April 2007
32. Scientific Committee on Emerging and Newly-Identified Health Risks. European Commission - Health & Consumer Protection DG - Directorate C: Public: Health and risk assessment opinion on the safety of medical devices containing DEHP plasticized PVC or others plasticizers on neonates and others groups at risk. http://ec.europa.eu/health/ph_risk/committees/04_scenihp/docs/scenihp_o_014.pdf
33. Yip E, Cacioli P. The manufacture of gloves from natural rubber latex. *J Allergy Clin Immunol* 2002;110(2 Suppl):S3-S14
34. Berekaa MM, Linos A, Reichelt R, Keller U, Steinbuchel A. Effect of pretreatment of rubber material on its biodegradability by various rubber degrading bacteria. *FEMS Microbiol Lett* 2000;184:199-206
35. Rahaman WA. Natural rubber as a green commodity. *Rubber Dev* 1994;47:13-16